

Mirsky's Fixative™



APPLICATION

Fixing Tissue Samples 105

- Saponin/Glyoxylate Fixative
- Preserves Immunohistological Activity
- Contains No Formaldehyde
- Fast-Acting
- No Toxic Buffers

Mirsky's Fixative is a superior fixing agent for use in immunohistological and immunocytological staining protocols. Mirsky's Fixative does not contain formaldehyde or glutaraldehyde. Therefore, it has considerably reduced toxicity and virtually no odor. Mirsky's Fixative is neutral, buffered, and isotonic (308 milli-osmol). Additionally, Mirsky's Fixative does not contain toxic or hazardous buffers such as cacodylate or barbital, and can be safely used for a wide variety of tissue samples. With Mirsky's Fixative, hardening or shrinkage of tissue is considerably reduced.

Double reactive sites afford excellent crosslinking properties while maintaining sample enzyme activity. Therefore, samples processed in Mirsky's Fixative for light microscopy can subsequently be used in electron microscopy procedures. Special buffer systems may be used in place of the buffer provided.

Method: Use as a replacement for formalin and/or glutaraldehyde fixatives in immunohistological and immunocytological staining protocols. Tissue size is unlimited as long as the sample is no thicker than 0.5 cm in at least one plane to assure uniformity of tissue penetration. Mirsky's Fixative is crosslinking and will therefore retain cellular morphology better than non-crosslinking fixatives such as formaldehyde. It is expected that the gross visual appearance will be different than in formaldehyde and tissue may seem "raw". This is due to double site binding of the fixative with a resultant reduction in tissue brittleness and shrinkage. Microscopic examination of tissue morphology will be noticeably improved. This material is intended to maintain enzyme and antibody activity. For best results, tissue sections of high digestive enzyme content should be thoroughly rinsed in saline solution before fixation (e.g. trypsin in intestinal samples). Once fixed, tissue may be retained in Mirsky's Fixative and ethanol and is stable indefinitely.

Mirsky's Fixative is normally distributed as a concentrated two bottle system, although high through-put laboratories may be interested in the single bottle ready-to-use format. The two bottle system is comprised of Mirsky's Fixative 10X Concentrate and Mirsky's Fixative 10X Buffer. To reconstitute to working strength, add 1 part Mirsky's Fixative 10X Buffer to 8 parts distilled or deionized water, mix thoroughly, then add 1 part Mirsky's Fixative 10X Concentrate and mix again. The material is now ready for use.

The two bottle system has a shelf life of 12 months while the ready-to-use format has a shelf life of 30 days.

Product Name	Cat. No.	Size
Mirsky's Fixative	HS-102	200 ml system 2 liter system
Mirsky's Fixative (ready-to-use)	HS-101	4X1 gal. 1X5 gal.