

Tissue Preparation

Whether you need to decalcify bone or to fix tissue to preserve fine structure, National Diagnostics has a product that will help. Whenever possible our products are designed to reduce hazards to the operator without compromising performance.



Calci-Clear™ Calci-Clear Rapid



Calci-Clear may be used to decalcify any tissue type. This solution of efficient chelating and sequestering agents removes and binds all calcium and other metal ions under treatment. Interference with subsequent staining is reduced while tissue distortion and structure loss is minimized.

Calci-Clear *Rapid* is a high-speed decalcifier which permits the processing of tissue samples in considerably less time than with other decalcifiers. Tissue distortion and structure loss are minimized while excellent staining characteristics are maintained.

Product Name	Cat. No.	Size
Calci-Clear	HS-104	1 quart
		1 gal.
		5 gal.
Calci-Clear <i>Rapid</i>	HS-105	1 quart
		1 gal.
		5 gal.

APPLICATION

Decalcifying Tissue Samples 107